

R.K. VIGYAN MAHAVIDHYALAYA
(Kalwar Jaipur)

NISHA SHARMA.

M.Sc.(Previous) Botany.

Roll no.-1259516.

$$R_f = \frac{\text{Distance from origin travelled by compound}}{\text{Distance of solvent front from origin}}$$

$$\text{Carotene} = \frac{11}{11.5} \text{ cm} = 0.96 \text{ cm.}$$

$$\text{Xanthophyll} = \frac{10.7}{11.5} \text{ cm} = 0.93 \text{ cm.}$$

$$\text{Chlorophyll 'a'} = \frac{3.5}{11.5} \text{ cm} = 0.30 \text{ cm.}$$

$$\text{Chlorophyll 'b'} = \frac{7.8}{11.5} \text{ cm} = 0.68 \text{ cm.}$$

1.

EXPT. NO.
DATE :

AIM : Separation of chloroplast pigment by paper chromatography technique.

MATERIALS : Spinach leaves, chromatography paper (Whatman No.1), mortar, pestle, separatory funnel, capillary tube, pencil, drier, benzene, distilled water, acetone.

PROCEDURE : The following are major steps.

- A). Preparation of paper :- 1. Cut chromatography paper (Whatman no.1) into square sheets to a size which would fit in easily in the chromatography jar.
2. Draw a pencil line $\frac{1}{2}$ inch above the bottom of each sheet.
- B). Preparation of extract :- Add 50ml of pre-cooled acetone to leaf pulp & stir well then filter through a fine linen cloth.
- C). Application of Paper chromatography :- 1. Put a spot of the extract with the help of a capillary tube $\frac{1}{2}$ inch from the left hand margin on the pencil line.
2. Allow the pigment spot to dry.
3. Fix the strip with clips to cork of tube or roll it around a glass at the top of the jar, so that its bottom just touches the solvent.
- D). Development :- 1. Benzene \rightarrow Acetone, 85:15
2. Pour sufficient solvent into chromatography jar filling about an inch from the bottom.

3. Place the spotted paper vertically, so that spot is just above vertical level.
4. Cover the jar and close the lid tightly.
5. Allow 1-2 hours for development.
6. The solvent shall move up the paper.
7. Remove the chromatogram when the solvent reaches the top and allow it to dry.
8. Since the colours of the pigments have examined the papers.


E). Result :- The sequence of pigments from top to the bottom shall be-

Carotene	- orange yellow.
Xanthophyll	- one or more yellow bands.
Chlorophyll 'a'	- blue green.
Chlorophyll 'b'	- yellow green.


3. Determination of the R_f values by following formula -

$$R_f = \frac{\text{distance from origin travelled by compound}}{\text{distance of solvent front to origin}}$$

Identify the compounds/components by comparing results with standardised values.



Demonstration of Paper chromatography.



Thin layer chromatography

EXPERIMENT—

2.

EXPT. NO.
DATE :

3.

AIM :- To isolate and identify different amino acids by paper chromatography methods.

Requirements :-

- i) Materials :- Known amino acid glycine, aspartic acid.
- ii) Chemicals :- Butanol, glacial acetic acid, ninhydrin, distilled water.
- iii) Glass wear :- Chromatography jar with lid, capillary tube.
- iv) Miscellaneous :- chromatography paper Whatman no.1, paper clips, iron stand, thread grease, pencil, scale.

Principle :- Chromatography is an analytical technique. It is based on two different phases.

1. Solid-Liquid position.
2. Liquid-Liquid position.

1. Solid-Liquid position :- It is based on difference in solute position between a stationary phase which is paper and their dissolution in mobile organic phase.

2. Liquid-Liquid position :- It is based on difference in solute position b/w a stationary aqueous phase which is tightly bound to cellulose fibres of filter paper and an organic mobile phase which is moving through paper by capillary action.

Different amino acids will be separated on the chromatography

paper depending on their affinity for the solvent, the mobile or stationary phase amino acid are color less. therefore they can't be detected with naked eyes for this detection a reagent called nin-hydrin is used.

It formed purple coloured complex called Ruhemann purple with which give yellow coloured complex.

Procedure :-

1. A chromatography jar with a lid was taken and chromatography paper is cut according to size of jar, such that it did not touch the sides of jar.
2. A line is then marked on the paper at a distance of 1.5" from the base.
3. Different points were marked, liquid & distantly on that line for different amino acids.
4. Different amino acids are loaded with the help of capillaries, taking over that capillary tubes are not mixed.
5. Loading is done 8-10 times after spot is dried.
6. Mixture of butanol, glacial, acetic acid & water made in ratio of 5:1:4 mixed thoroughly and poured into jar.
7. The mixture is kept in the jar for 5-10 hours min.



8. Then the chromatography paper is folded in the form of a cylinder and kept in the jar.
9. The lid was placed after applying grease and the solvent was placed after it & allowed to run undisturably for one hour.
10. After one hour paper is removed and solvent level is marked.
11. Paper is dried and 1% ninhydrin prepared in acetone is separated with the help of atomisers.
12. Paper is then kept in oven for colour development (80°C).
13. Paper is taken out different amino-acids marked.

Result :- Given amino acid is aspartic acid.

- Precaution :-
1. Solvent must be prepared accurately.
 2. The capillaries of different amino acids should not be mixed.
 3. Apparatus should be air tight.
 4. The chromatography paper should not touch the water of jar.

AIM :- To extract and chemically separate chloroplast pigment by liquid-liquid chromatography.

Requirements :-

- (i). Materials :- Fresh leaf extract (spinach).
- (ii). Chemicals :- 80% acetone, Petrolatum ether, 92% methyl alcohol, dimethyl ether, 30% Methemolic KOH, distilled water.
- (iii). Glassware :- Separatory funnel, Bucher's funnel for filtration, test tube, beaker & funnel.
- (iv). Miscellaneous :- Filter paper, Pesto and marker.

Principle :- The basis of chemicals separation of chloroplast pigment is the relative immisibility of various solvent and the differential solubility of solvents in different solvents.

Both these properties depends upon the polarity of the solution and solvent molecules. A polar solute is more readily soluble in a polar solvent than in a non-polar solvent while the reverse is true for non-polar solute.

Thus of four major pigments of chloroplast chl-a, chl-b, carotene and xanthophyll, Xanthophyll being relatively more polar than chl-a and carotenoids either captures chl-a and carotenoids.

Procedure :- 1. Fresh leaves of spinach are ground with

S.no.	Concentration (mg)	optical density.
1.	.2 X 200 mg =	40
2.	.4 X 200 mg =	80
3.	.6 X 200 mg =	120
4.	.8 X 200 mg =	160
5.	.10 X 200 mg =	100

EXPERIMENT

EXPT. NO.
DATE:

80% acetone.

2. When it was deep green coloured filter it through filter paper.
3. filter extract again through Buchber's funnel.
4. The filter contains chloroplasts pigments.
5. In a separatory funnel, Put 50ml of Petroleum ether.
6. Gently coated the funnel till upper layer was deep green. The upper layer is petroleum ether layer and lower is acetone layer.
7. Acetone extract is added to it.
8. Distilled water is added to the funnel along the sides again until the upper layer is quite green and lower layer is colourless.
9. Gradually collect the lower acetone layer in a beaker and discard it.
10. To the clean petroleum ether layer added 50 ml of 92% methyl alcohol.
11. Mixed thoroughly by rotating gently till two layers were obtained

in the separating funnel.

12. The upper petroleum ether methyl alcoholic layer in to one beaker and upper petroleum ether layer in to another beaker.
13. Took the methyl alcoholic beaker in a separating funnel and add 50 ml of diethyl ether and mixed it by rotating.
14. Now add distilled water taking small at a time down the slide of the funnel till two layers of liquid appear.
15. The upper one is diethyl ether layer and lower is methyl alcoholic layer.
16. Discard the lower methyl alcoholic layer should be colourless.
17. Put some of diethyl ether layer in a building tube and carefully added 15 ml of fresh prepared 30% methanolic KOH.
18. The lower layer of methanol is formed having chl - b and upper diethyl ether layer has Xanthophyll.
19. To the petroleum ether layer which is poured into the another beaker is added 15 ml of 30% methanolic KOH.
20. Two layers are formed lower methanolic layer of chl a and upper petroleum ether layer of carotenoids.



EXPERIMENT

EXPT. NO.

DATE :

9.

Result :- four chloroplast Pigments.

Chl 'a' - blue green.

Chl 'b' - yellowish green.

Xanthophyll - dark yellow.

Carotene - orange yellow.

Precaution :- 1. Grind the leaves well.

2. The separating funnel has to be rotated gently.

3. Avoid the funnel of methyl alcohol.

AIM :- To prepare standard curve of soluble protein and calculate the quantity of Protein in given sample.

Requirements :- 1). Materials:- Beaker, test tube, stand, calorimeter.

Reagent A - 20% sodium carbonate (Na_2CO_3) + 0.1% sodium hydroxide (NaOH).

Reagent B - 0.5% copper sulphate + 1.0% potassium sodium tartate.

Reagent C - 50ml of reagent A + 1ml of reagent B prior to their use.

Reagent D - ~~foline - icatcaes~~ reagent protein solution (stock standard).

Procedure :- 1. Procedure for making standard curve working standard \rightarrow 10ml of the stock solution was diverted to 50ml with distilled water. A series of volume from 0.1 to 1.0 ml of this standard 30ml gives a concentration range from 120 mg - 200 mg protein ^2 was proceeded as that sample.

2. Procedure for estimation of Protein in unknown sample.

a). Extraction :- for extraction of Protein from the sample under ~~experimentation~~, 500 mg of sample is macerated.

Principle \rightarrow In the acidic environment of the reagent, Protein binds to the comassie dye. This results in a spectral shift from reddish brown form of the dye absorbance Max. at 465nm to blue form of the dye. Development of colour in coomassie dye

based protein assays has been associated with the presence of certain basic amino acids in the protein. Vander walls forces and hydrophobic interactions also participate in the binding of the dye by protein, the no. of coomassie dye ligands bounds each protein molecule is approximately proportional to the number of positive charges found on the protein free amino acids, Peptide and low molecular weight coomassie dye reagent. In general the mass of a peptide on protein must be at least 3,000 daltons to be assayed with this reagent. In some applications this can be an advantage. The residue was re-centrifuged and washed with 1M cold trichloroacetic acid centrifugation was performed again and the protein fraction was dissolved in 1M NaOH. This solution was made to a known volume and kept for 1 hour incubation at room temp.

(b). Estimation :- To 0.1 ml of dissolved residue 5 ml of reagent was added and allowed to stand for 10 min. To this 0.5 ml of folin ciocotteae reagent is rapidly added mixed thoroughly and incubated at room temp for 30 min. a blue colour developed. The O.D was measured at 660 nm. The amount of Protein is calculated with standard curve prepared from bovine serum albumin, Protein control is expressed in terms of $\mu\text{g}/\text{mg}$ fresh weight of the tissue. The result presented are average of three replicate.

Result :- The concentration of unknown protein sample is measured at 660 nm is 40 mg/ml.

Precaution → Glassware should be cleaned properly.

AIM:- Measurement of chlorophyll a and chlorophyll b in C₃ and C₄ plants by spectrophotometer.

Requirement :-

- (a). Material :- Grass leaves (cynodon - C₄) (spinach and bougainvillea - C₃ plants).
- (b). Chemicals :- Acetone 80%.
- (c). Glass wares :- Beaker, measuring cylinder, test tube.
- (d). Miscellaneous :- Spectrophotometer.

Principle :- C₃ and C₄ plants are group of plants named due to respective photosynthesis is operating in them. These two groups of plants differ in rate of photosynthesis and also in their adaptability of towards particular environmental conditions.

Rate of Photosynthesis is less in case of C₃ Plants of they consists of 2 types of chloroplasts granular and agranular types therefore rate of photosynthesis is twice as compared to C₄ plants because C₄ plants are capable to converting solar energy into dry matter because of all these reasons chlorophyll content in C₄ plants also the ratio of chlorophyll a and b is more in case of C₄.

$$C_4 \text{ plants} = a:b = 3:9:1$$

$$C_3 \text{ Plants} = a:b = 2:8:1$$

Absorption max. for chl b is 645 nm.

Procedure :- 1. Grind 1 gm of grass leaves (C₄) and Bougainvillea (C₃) leaves separately with small amount of 80% acetone.



Calculation for C₃ Plants :-

$$\text{For chl}'a' = (12.7 \times A_{663} - 2.69 \times A_{645}) \times \frac{V}{1000 \times \text{wt}} \text{ mg/g fresh wt.}$$

$$= (12.7 \times 6.7 - 2.69 \times 30) \times \frac{50}{1000}$$

$$= (8.609 - 80.7) \times \frac{50}{1000}$$

$$= 7.702 \times \frac{50}{1000}$$

$$= 0.389 \text{ mg/g fresh wt.}$$

$$\text{For chl}'b' = (22.9 \times A_{645} - 4.68 \times A_{663}) \times \frac{V}{1000 \times \text{wt}} \text{ mg/g.}$$

$$= (22.9 \times 30 - 4.68 \times 6.7) \times \frac{50}{1000}$$

$$= (6.87 - 3.13) \times \frac{50}{1000}$$

$$= 3.74 \times \frac{50}{1000} = 1.87 \text{ mg/g fresh wt.}$$

$$\text{Total chl} = (20.2 \times A_{645} + 8.03 \times A_{663}) \times \frac{V}{1000 \times \text{wt}} \text{ mg/g.}$$

$$= (20.2 \times 30 + 8.03 \times 6) \times \frac{50}{1000 \times \text{wt}}$$

$$= (6.06 \times 5.37) \times \frac{50}{1000}$$

$$= 11.4337 \times \frac{50}{1000}$$

$$= 0.57167 \text{ mg/g fresh wt.}$$

EXPERIMENT

EXPT. NO.
DATE :

2. filter the above sample.

3. Take optical density (O.D) of these and sample at their absorption maximum separately i.e. at 663 nm and 645 nm and calculate chlorophyll content by the following formulae.

$$\text{Chl}'a' = 12.7 \times A_{663} - 2.69 \times A_{645} \times \frac{V}{1000 \times \text{wt(gm)}} \text{ mg/g.}$$

$$\text{Chl}'b' = 22.9 \times A_{645} - 4.68 \times A_{663} \times \frac{V}{1000 \times \text{wt(gm)}} \text{ mg/g.}$$

$$\text{Total chlorophyll} = 20.2 \times A_{645} + 8.03 \times A_{663} \times \frac{V}{1000 \times \text{wt(gm)}} \text{ mg/g.}$$

where $A_{663} = 0.0$ at 663 nm.

$A_{645} = 0.0$ at 645 nm.

V = Volume of extract in ml.

wt = fresh wt. of the sample.

Result :- Chlorophyll content in higher in C₃ plants than C₄ plants.

In C₃ Plant - chl}'a' = 0.4515 mg/g fresh wt.

chl}'b' = 0.2937 mg/g fresh wt.

total chl = 0.745 mg/g fresh wt.

In C₄ Plants \rightarrow Chl'a' = 0.3851 mg/gm fresh wt.

Chl'b' = 0.187 mg/gm fresh wt.

total chl = 0.5716 mg/gm fresh wt.

Precaution :- 1. Grinding should be uniform.

2. Spectrophotometer should be turned on before 15 min
the experiment.

AIM :- To study the activity of Enzyme Catalase.


Requirement :- Potato tubers, knife, Petridishes, test tubes, test tube holder, test tube stand, spirit lamp, hydrogen Peroxide.

- Procedure :-
1. cut small and thin pieces of potato tubes.
 2. Place the pieces in petridish filled with water.
 3. Take a few small pieces and put them in two separate test tubes marked A and B fill both the test tube with appropriate amount of water.
 4. Keep the test tube A in the stand boil the piece of Potato tuber placed in tube B.
 5. Drain out water from both the test tube, now add hydrogen peroxide solution so as to completely immerse the pieces of Potato tubers.
 6. Note the changes occurring in both the test tubes.

Result :-

1. Evolution of bubbles is seen in test tube A.
2. No bubbles seen in test tube B.

Conclusion :- Catalase brings about decomposition of H_2O_2 into water and oxygen.



- O_2 evolution, therefore is an O_2 reaction of activity of catalase.
- The oxygen bubbles are not evolved in test tube B. This is because enzyme is destroyed when potato slices are boiled Hence, there is no enzymatic activity.

AIM :- To study the phenomenon of plasmolysis.

Requirements :- Onion, Safety blade, Sugar solution, cover slips, water, slides, microscope etc.

Procedure :- 1. Peel off a small segment from the onion surface. This can be done by tearing the surface obliquely with a single jerk of scraping it with safety blade.

2. Mount the bell in a drop of water on a slide and then place a cover slip observe under the microscope. Draw the protoplasm and set the preparation it is called A. Take another bell, cut and similarly mount pieces in a drop of sugar solution observe preparation under the microscope Draw the boundaries of Protoplasm, let the preparation be called as B.

Results :- 1. In condition A, the cell structure can be seen clearly. The cells are turgid and protoplasm is closely pressed against cell wall.

2. In condition B, when a concentrated sugar solution is used. the cell contents move appreciably away from the cell wall, leaving a considerable space b/w the wall and sap.

Conclusions :- 1. Preparation A shows normal condition where cell sap presses the protoplasm against the cell wall which is slightly inflated. The cell wall is called turgid.

2. Preparation B, with the increase in the concentration

EXPERIMENT

EXPT. NO.

DATE :

17.

of sugar solution outside leaves from the cell comes out and the space b/w cell wall and the contents increases. finally due to continued exosmosis the cell contents shrink and collect on one side. such a cell is called as plasmolysed.

AIM :→ To perform qualitative test for protein.

Requirements :→ Test tubes, water bath, pea seeds, spirit lamp, conical flask, test tube holder, chemicals such as HNO_3 , H_2SO_4 , Ninhydrin, copper-sulphate, sodium hydroxide, ammonium hydroxide, millon's reagent A & B etc.

Procedure :→ Take pea seeds and mash in mortar with 80% alcohol. This residue is discarded and filtrate is used for performing various test in test tubes.

1. Hellig's test :→ In this test yellow colour occurs due to conversion of protein to meta protein.

Experiment	Observation	Inference.
1. filtrate + few drops of HNO_3 + boil + add few drops of lead acetate.	white colour turns to yellow ppt.	Protein is present.

2. Hopkin's and test :→ In the presence of acid many aldehydes condense with indole ring of tryptophane to yield coloured product. Presence of tryptophane is indicated by the appearance of violet ring at interphase.

Experiment

Observation

Inference.

Protein extract +
1 ml of conc H₂SO₄
+ 2 ml of glacial
acetic acid.

Reddish brown
colour appears.

Protein Present.

3. Nin-hydriin test :→ Nin-hydriin solution reacts to reduce Nin-hydriin
reacts to produce purple blue compound which is stable in
nature.

Experiment

Observation

Inference.

Protein extract +
nin-hydriin soln +
warm.

Blue colour appear

Protein present.

4. Xanthoprotein test :→ This test is performed for aromatic and heterocyclic amino acid presence in the protein amino acid which contain aromatic nucleolus forms yellow nitro derivative on heating with conc HNO₃ when cooled and 20 × NaOH in sufficient amount is added then yellow colour changes to bright orange in alkaline conditions give positive result.



Experiment	Observation	Inference.
Protein extract + conc HNO ₃ + heat + few drops of 20% conc ammonium hydroxide.	white ppt is formed yellow colour developed and it changes to orange.	may be Protein Present Protein Present.

5. Bi-urete test :- Alkaline cu soy reacts with components containing two or more proteins bonds to give violet colour complex. The depth of the colour is a measure num of peptide bonds present in protein.

Experiment	Observation	Inference.
Protein extract 1 ml. NaOH (20%) + few drops of cu soy + Now add millen's reagent (A+B) & Boil.	Violet colour develop Reddish colour appear	Protein may be Present. Protein present.

Conclusion :- The change in colour indicates the presence of Protein.

AIM : → To perform the qualitative test for organic acid.

Requirements : → Paste and mortar, test tube, test tube holder, Spirit lamp, many chemicals lemon, tamarind, tomato.

Procedure : → 1. Test for tataric acid : → Grind the tamarind in distilled water filter and neutralize with Ammonium hydroxide (NaOH).

Experiment	Observation	Inference.
A. Extract + add 5% CaCl_2 .	PPT is formed	tartaric acid present.
B. Extract + few drops of Ammonical AgNO_3 (0.1N) + 8ml NH_4OH + water + warm.	Silver mirror is formed at the bottom of test tube.	tartaric acid present.

2. Test for malic acid : → Take tomato juice and add water and neutralize with NH_4OH .

Experiment	Observation	Inference.
A. Extract + 5% CaCl_2 + boil and cool + equal amount of ethanol (95%).	white PPT is formed.	malic acid present.

B. Extract + few drops of lead acetate (5%). Now add warm acetic acid.	white ppt appears Ppt is soluble	Malic acid present.
C. Extract + ferric chloride FeCl_3 .	yellow colour is obtained	Malic acid confirmed.

3. Test for citric acid :→ Take lemon juice and add little amount of distilled water neutralize with NH_4OH ammonium hydroxide.

Experiment	Observation	Inference.
Extract + 5% CaCl_2 + heat + acetic acid.	No. ppt formed white ppt appears.	citric acid present.

4. Test for ascorbic acid :→

Experiment	Observation	Inference
lemon juice dropped on filter paper and treated with KMnO_4 .	filter paper is decolorised due to oxidative reaction of ascorbic acid.	Ascorbic acid present.

Conclusion :→ The change in colour indicate the presence of organic acids.

AIM :→ To demonstrate the presence of Anthocyanin and tannin in the given sample (Plant material).

Requirements :→ Test tubes, Paste and Morter, bouganvillia, tea leaves, NH_4OH , potassium ferrocyanide, HNO_3 , FeCl_3 , NaOH , HCl , KMnO_4 etc.

Procedure :→ 1. Test of tannin :→ Take 5gm of tea leaves and boil in water for 30 minutes. Now filter and extract is taken into three test tubes.

Experiment	Observation	Inference.
1. Extract + NH_4OH + Potassium ferrocyanide	Red colour appears.	Tannin is Present.
2. Extract + 1 drop of HNO_3 .	Yellow colour appears.	Tannin is present.
3. Extract + Few drops of FeCl_3 .	Blue - colour formed.	Tannin is present.

2. Test for anthocyanin :→ Macerate bouganvillia flowers in mortar with 1× HCl in ethanol. Take the sample by filtration in test tubes (three).

EXPERIMENT

EXPT. NO.
DATE :

24.

Experiment	Observation	Inference.
1. Extract + dil HCl.	Red colour appears.	Anthocyanin present.
2. Extract + 1% al- coholic Fe(II)	Yellow colour appears.	Anthocyanin present.
3. Extract + Ag-NaOH.	Blue colour formed.	Anthocyanin present.

Conclusion :- The colour formed indicates the presence of tannin and anthocyanin.

AIM :→ To study the permeability of plasma-membranes using Organic Solvents.

Requirements :→ Beet roots, distilled water, test tubes, test tube holder, test tube stand, cork, beaker, ethanol, benzene, Petroleum etc.

- Procedure :→
1. cut slices or cylinders of beet root by razor or cork borer.
 2. Repeatedly wash the slice with distilled water until the pigment stop diffusing into water.
 3. fill the test tubes A with distilled water, B - with ethanol, C with Benzene, D with Petroleum ether.
 4. Place a few slices of beet root or cylinders in each test tubes.
 5. Allow the test tubes to remain as such for about 30 min.

Result :→ Since the red anthocyanin pigment diffuses in a order ethanol, benzene, water and Petroleum ether, Permeability is lowest in petroleum ether and highest in ethanol.

Conclusion :→ The appearance of red color in the test tubes containing ethanol is due to diffusion of anthocyanin pigment from inside the cell to extract medium in a living plant cell selectively.

Permeable cell membrane prevents out-ward diffusion of anthocyanin. Toxic substances like ethanol kill the cell therefore the cell membrane also which

EXPERIMENT

EXPT. NO.

DATE :

26.

Loses its selective permeability and allows anthocyanin to diffuse out.

Precaution :- 1. Root should be mashed thoroughly.
2. Cylinder should be of uniform size.
3. Time should be constant for all test tubes.

AIM :→ To study the effect of temperature on permeability of plasma membrane.

Requirements :→ Beet roots, cork borer, distilled water, test tube, test tube holder, test tube stand, tripod stand, beaker heating plate etc.

Principle :→ It is based on the Principle that red colour pigment anthocyanin is present in beet root. Anthocyanin pigment which are water soluble. In normal condition anthocyanin doesn't come out of the membrane. Temperature affects the permeability of the membrane and in turn affects the concentration of anthocyanin that come out of the membrane.

- Procedure :→
1. Cut slices of beet root using cork borer.
 2. Wash the slices repeatedly with distilled water till pigment stops diffusing into the water.
 3. Fill the test tube (A) with tap water, (b) with cold water, (c) with hot water (d) with boil water.
 4. Place a few slices of beet root in each test tube.
 5. Allow the test tube to remain as such for about 30 min.

Result :→ The red pigment anthocyanin diffuses in order boil water, hot water, tap water then cold water.

→ Permeability is lowest in cold water and highest in boiled water.

EXPERIMENT

EXPT. NO.
DATE :

27.

Conclusion :→ The diffusion of anthocyanin is the result of death of cell membrane due to high temperature.

Precaution :→ 1. Root slices should be washed separately.
2. Slices should be of uniform size.
3. Time should be constant for all test tubes.

Reagent	Gram (-)	Gram (+)
None heat fixed cell	◦ //	◦ ◦ //
Crystal violet (30 sec)	◦ //	◦ //
Gram iodine (20 sec)	◦ //	◦ //
Ethyl alcohol (20 sec)	◦ //	◦ //
Sulphur (20 sec).	◦ //	◦ //

EXPERIMENT

13.

EXPT. NO.

DATE:

24.

AIM : → Gram staining bacteria.

Requirement : → 24 hour old culture of bacteria, staining tray, inoculating loop, slide blotting paper, spirit lamp, microscope.

Reagent : → Crystal violet, gram iodine soln, ethanal, Sabazine.

Theory : → Gram staining was developed by Dr. Henery Christian gram, it's useful staining method for identify and classify bacteria into two major group-

1. Gram + 2. Gram -

- Crystal violet → Primary stain.
- Iodine solution → mordants
- Alcohol → decolorizing agent.
- Sabazine → counter stain.

The bacteria that retain to primary stain are called gram +, which loose it called gram - bacteria.

Principle : → The different stain responses to gram staining can be related to chemical and physical differentiate in their cell wall.

The gram + cell wall is this composition mouth layer higher liquid content addition to protein and

EXPERIMENT

EXPT. NO.
DATE :

30.

mucopeptide. The hidden layer amount of lipid is being readily dissolution by alcohol leading to be formation a large pores in the cell wall.

In counter the gram \ominus cell wall are thick and chemically simple.

make a thin smear air dry and heat the smear to thin it.

AIM :→ To detect the presence of starch.

Requirements :→ Starch, Iodine solution, hydrochloric acid, benedict's solution sodium hydroxide or sodium-carbonate (NaOH or Na_2CO_3), test tube, test tube holder, Spirit lamp etc.


Procedure :→ Following two methods are used to detect the presence of starch.

Method I → i). Take a pinch of starch in a test tube and add water.
ii). Add few drops of Iodine solution to it.

Method II → i). Take a pinch of starch in a test tube and add water to it.
ii). Add hydrochloric acid into it and boil the solution.
iii). Make this solution alkaline by adding NaOH or Na_2CO_3 .
iv). Now add Benedict's solution and observe.

Results & conclusion :→ i). In first experiment blue black colour appears.

ii). Red precipitate is formed with Benedict's solution. The test indicate the presence of starch.



(Apical meristem culture) (virus free plant culture).

EXPERIMENT

15.

EXPT. NO.
DATE:

32.

Comments \Rightarrow Virus free culture.

Manal and Merlin (1952) developed the thin of meristem culture for in vitro virus enculturation of Dahlia.

Colages marel (1965) was the pioneer in applying shoot tip culture for micro propagation of orchid Typhlidium.

When diseased free plants are the objective meristem tips should be cultured.

Meristem tips b/w 0.2 and 0.5 mm most frequently produce virus free plants and this method is referred to as meristem tip culture.

The apical meristem is usually a dense a dense tissue located at the extreme tip of a shoot and measures ~ 0.1 mm in diameter and ~ 0.25 to 0.3 mm in length.

Meristem and shoot tip is cut or isolated from stem by applying a V-shaped cut with a sterile knife.

Requirement ie flask / culture tubes etc.

- \Rightarrow fold the act piece from the centre
- \Rightarrow Roll the sheet from any of ends opposite to the folded side has cotton fibres that energy out freely.

- ⇒ Insert the folded and rolled side onian the neck / mouth of a flask culture tube and the other end of plug bears cotton fibres.
- ⇒ If it is tight open the plug and remove same cotton from one side and read just or earlier.
2. Cut rolled and fold plug ⇒
Cut appropriate place of cotton sheet with the help of scissors.
- ⇒ Roll the piece from one end.
- ⇒ Roll the rolled piece from the centre.
- ⇒ Insert the folded one in the neck mouth of flask culture tube and the other end bearing cotton fibres out in air.

Result :- The isolated colony of desired microbes cut the side of last streak on these plates will be observed.

WJL
18/10/2018