

NAME:- MUKESH MEEL

CLASS:- M.Sc. (Previous)

SUBJECT:- PHYSICAL CHEMISTRY

COLLEGE:- R. K. VIGYAN (PG) COLLEGE

KALVAR, JAIPUR

ROLL.No:- 1285828

Sl No.	Name of the Experiment	Page No.	Date of Experiment	Date of Submission	
(1)	determine the solubility and solubility product of sparingly soluble salt $BaSO_4$		Silvany		
(2)	Determine the strength of given strong acid conductometrically using standard strong base solution.		↑		
(3)	Determine the strength of a given strong acid (HCl) conductometrically using standard base NH_4OH				
(4)	Determine the strength of given weak acid g conductometrically using strong base $NaOH$.				
(5)	calculate the rate constant for the reaction of ethyl acetate			Silvany	
(6)	determine the triangular phase diagram of CH_3COOH , HCl , H_2O				

object:- Determine the solubility and solubility product of sparingly soluble salt BaSO_4 .

Principle:-

Conductance measurement may be employed in determination of solubility of a sparingly soluble salt like barium sulphate. Since the saturated solution of a sparingly soluble salt is very dilute $c_{\text{BaSO}_4} < 10^{-3} \text{ N}$ so it may be regarded as completely ionised in soluble state its equivalent conductance is that infinite dilution obtained by assumption of ionic conductance of the constituent ions.

In the present case

$$\lambda_{\infty} = \lambda_{\infty}^+ + \lambda_{\infty}^-$$

As the specific conductance of water must be taken in to account.

Teacher's Signature : _____

Calculation:-

Correct the specific conduction of the solution of a sparingly soluble salt by subtracting the specific conductance of water.

$$K = K_{\text{solution}} - K_{\text{water}}$$

$$= 0.375 - 0.372$$

$$= 0.003 \times 10^{-3} \text{ S}$$

The solubility of a sparingly soluble salt is.

$$S = \frac{1000 \times K}{\lambda_{\infty}} \times \text{equivalent weight}$$


$$= \frac{0.003 \times 10^{-3} \times 1000}{146.16}$$

$$= 0.0024 \text{ gm/liter}$$

solubility product of $\text{BaSO}_4 = S^2$

$$= 0.0024 \times 0.0024$$

$$= 0.00000484$$



$$K = K_{\text{solution}} - K_{\text{water}}$$

$$\rho = K_{\text{solution}}$$

$$\frac{K \times 1000 \text{ gm equivalent/liter}}{\lambda_{\infty}}$$

procedure:- (1.) Determine the cell constant of conductivity cell at 25°C.

(1.) measures the conductance of water used in preparation of solution.

2.) prepare fresh barium sulphate solution with hydroxide acid wash the precipitate thoroughly with conductivity water for several times to remove soluble and adhered impurities present if any.

place the solution well washed salt in a thermostat at 30°C stir it well after half an hour pipette out the clear solution in the beaker and measured the conductance are obtained.

Teacher's Signature : _____

observation:-

Temperature = 30°C

S.N.	observed conductance of water (milliS)	observed conductance of BaSO ₄ (milliS)
1.	0.371	0.374
2.	0.375	0.370
3.	0.370	0.381
	Average = 0.372	0.375

Result:-

solubility of the sparingly soluble salt BaSO₄ is 0.0024 gm/liter and solubility product is 0.00000484.

Surendra
23/4/18


Teacher's Signature : _____

object:-

Determine the strength of given strong acid conductometrically using standard strong base solution.

principle:-

Titration of strong acid like hydrochloric acid with a strong base like sodium hydroxide initially the value due to chloride and hydrogen ions as NaOH solution is gradually added to it. The hydrogen ions are replaced by ions thus decreased.



It goes on decreasing until HCl is completely neutralised i.e. the equivalence point reached. Now after this point any subsequent addition of sodium hydroxide.

solution result in increase of conductance due to Na^+ ions and OH^- ions. A plot of conductance vs volume of titrant gives two straight lines the point of intersection of extrapolated lines will be the end point of titration and gives the value of alkali required for the neutralisation of the acid.

Requirement:-

conductivity meter conductivity cell
microburets the pipette brokes.

chemicals:- $N/10$ NaOH and unknown till

procedure:- (i) $N/10$ NaOH standard solution

is prepared and filled in microburette.

(ii) 20ml given hydrochloride acid is taken in
baker (100ml)

(iii) set the conductometer at cell constant and top
conductance of given solution of till.

(iv) Now 0.5ml NaOH added from microburette and conductance
is measured this process is repeated after each success
in addition of alkali.

(v) Now graph is plotted and point found it

observation:-

strength of NaOH = $N/10$

Volume of acid = 20ml

Calculation:-

Strong acid and strong base.

$$N_1 V_1 = N_2 V_2$$

$$N_1 = \frac{N_2 V_2}{V_1}$$

N_1 = Normality of acid

V_1 = Volume of acid

N_2 = Normality of base = $\frac{1}{1} = 1$

V_2 = Volume of base = 1.8

$$N_1 = \frac{1}{10} \times \frac{1.8}{20} N$$

$$= 0.009 N$$


Concentration of given hydrochloride acid

$$= N_1 \times \text{equivalent weight}$$

$$= 0.009 \times 36.5 \text{ g/L}$$

$$= 0.3285 \text{ g/L}$$

1 cm³ 1 big black 0.2 ml
2 cm³ 1 big black 0.4 ml
mls



S.No.	Volume of alkalic NaOH added (in ml)	Conductance of the solution.
1.	0.0	1.452
2.	0.2	1.321
3.	0.4	1.223
4.	0.8	1.084
5.	1.0	0.998
6.	1.2	0.592
7.	1.4	0.523
8.	1.6	0.385
9.	1.8	0.486
10.	2.0	0.608
11.	2.2	0.662
12.	2.4	0.753
13.	2.6	0.825
14.	2.8	0.926
15.	3.0	1.002
16.	3.2	1.099
17.	3.4	1.86
18.	3.6	1.239
19.	3.8	1.357
20.	4.0	1.442
21.	4.2	1.573
22.	4.4	1.697
23.	4.6	1.800
24.	4.8	1.756
25.	5.0	1.827

Result:- The strenght of given hydrochloric acid is
0.3285 gm/litter.

Teacher's Signature : _____


object:-

Determine the strength given strong acid (HCl) conductometrically using standard base (NH₄OH).

principle:-

When a titration is carried out between strong acid and weak base like HCl and NH₄OH respectively the conductance decreases due to the replacement of H⁺ ions by slow moving NH₄⁺ ions until the equivalence point is reached after the end point and further

addition of NH₄OH does not cause any appreciable change in conductivity. The graph becomes horizontal parallel to volume axis since the excess of equimolar ammonia is not appreciably ionised in the presence of ammonia chloride.



Requirements - conductivity meter microburette
pipette, beaker

Chemicals -

$\frac{N}{10}$ NH₄OH (standard solution)

and known HCl

Conclusion:-

Weak base and strong acid

$$N_1 V_1 = N_2 V_2$$

$$N_1 = \frac{N_2 V_2}{V_1}$$

N_1 = Normality of acid = ?

V_1 = Volume of acid = 20

N_2 = Normality of base = 0.1

V_2 = Volume of base = 1.8

$$N_1 = \frac{1}{10} \times \frac{1.8}{20} N$$


Concentration of hydrochloric acid

$$= N_1 \times \text{equivalent weight}$$

$$= 0.009 \times 36.5$$

$$= 0.3285 \text{ gm/liter}$$

radius. 1 big black = 0.2 ml
y axis = 1 big black = 0.1 ml
ml has



procedure:-

- (i) $N/10$ NH_4OH standard solution is prepared.
- (ii) 20ml HCl and is taken in the beaker.
- (iii) set the conductometer on cell constant.
- (iv) Take the conductance of cell.
- (v) Now add 0.2ml NH_4OH and conductance is measured addition of alkali.
- (vi) Now graph is plotted and end point is found out.

observation:-

Strength of NaOH = 0.1N
volume of acid = 20ml

Teacher's Signature : _____

S.No.	Volume alkali NH_4OH (added in mL)	
1.	0.0	0.688
2.	0.2	0.673
3.	0.4	0.652
4.	0.6	0.633
5.	0.8	0.614
6.	1.0	0.592
7.	2.0	0.576
8.	2.2	0.557
9.	2.4	0.537
10.	2.6	0.517
11.	2.8	0.517
12.	3.0	0.518
13.	3.2	0.517
14.	3.4	0.516
15.	3.6	0.519
16.	3.8	0.520
17.	4.8	0.517

Result:-

The strenght of given hydrochloric acid is
 0.3285 gm/litter .

Teacher's Signature : _____

object:-

Determine the strength of given weak acid and CH_3COOH and conductrometrically using standard strong base (NaOH).

Principle:-

If titration is carried out between a weak like acetic acid and strong base ex:- sodium hydroxide solution the conductance due to acid is very low on account of its poor dissociation. It first decrease using to the formation of weak acid but soon after the in a increase in conductance after complete neutralization of the acid the conductance \uparrow rapidly due to increase in OH^- ions.

**Requirement:-**

conductivity meter, microburette, pipette, beaker.

chemicals:-

$\frac{N}{10}$ NaOH (standard solution and

Unknown, CH_3COOH (acetic acid)

Teacher's Signature : _____

calculation :-

Strong acid and weak base

$$N_1 V_1 = N_2 V_2$$

$$N_1 = \frac{N_2 V_2}{V_1}$$

N_1 = Normality of acid = ?

V_1 = volume of acid = 20

N_2 = Normality of base = 1

V_2 = volume of base = 1.1

$$N_1 = \frac{1}{10} \times \frac{1.1}{2.0}$$
$$= 0.0055N$$


Concentration of given acetic acid

= $N \times$ equivalent water

$$= 0.0055 \times 60$$

$$= 0.33 \text{ gm/l}$$

x-axis = 4 small black = 1 ml
y-axis = 1 big black = 0.5



procedure:-

- (i) $N/10$ NaOH solution is prepared and filled in microburette.
- (ii) 20ml CH_3COOH acid is taken in beaker.
- (iii) set the cell constant of conductivity meter take the conductance of CH_3COOH .
- (iv) Now add 0.2ml NaOH and conductance is measured.
- (v) This process is repeated after each successive addition of alkali.

observation:- Normality of NaOH = $N/10$

Volume of acid = 20ml

S.N.	Volume alkali (NaOH) add (in ml)	conductance of solution (ml)
1.	0.0	0.642
2.	0.2	0.610
3.	0.4	0.516
4.	0.6	0.623
5.	0.8	0.752
6.	1.0	0.871
7.	1.2	0.990
8.	1.4	

Teacher's Signature : _____

9.	1.6	1.123
10.	1.8	1.267
11.	2.0	1.367
12.	2.2	1.491
13.	2.4	1.600
14.	2.6	1.738
15.	2.8	1.851
16.	3.0	1.951
17.	3.2	2.050

Result:-

The strenght of given acid is 0.33 gm/litter.

~~S. Ahmed~~
23/4/18

object:-

Calculate the rate constant for the reaction of saponification of ethyl acetate by sodium hydroxide solution conductometrically.

Principle:-

When during a reaction ions are produced and used of ~~are~~ replaced by some other ions then the conductance of the solution under goes change as the saponification of the acetate by sodium hydroxide the stoichiometry is represented as.



In the saponification reaction the fast moving of OH^- ions are replaced by slow moving CH_3COO^- ions. Initially the NaOH solution and finally it will be due to ions formed by an equivalent amount of sodium are produced as the reaction process the conductance of reaction mixture thus decreased if the reaction conductance of reaction

is carried out with same initial concentrations. The reactions then the rate constant can be calculated using integrated rate law for the second order reaction.

Teacher's Signature : _____

$$k = \frac{1}{t} \frac{n}{a(a-n)}$$

(k = second order rate constant)

Requirement:-

(a) conductivity meter conductivity cell stop watch volumetric beaker burette pipette etc.

(b) 0.05 NaOH phenolphthalein ethyl acetate - acetic acid.

Procedure:-

(i) 0.05 NaOH is prepared and standardized.

(ii) solution A:- 20ml 0.05 NaOH and make it up to 100ml now this is the 0.01N solution A.

(iii) solution B:-

(a) first 0.05N NaOH is taken in a beaker and add 50ml water now it is 10ml.

(b) 5ml ethyl acetate in other beaker and add 25ml water now in 30ml now both (a) and (b) solution to make solution 100ml co-value is taken by a solution with the help of conductometer.

Teacher's Signature : _____

Calculation:-

$$K_1 = \frac{1}{0.01} \times \frac{C_0 - C_t}{C_t - C_\infty} = \frac{1}{0.01} \times \frac{1.412 - 1.167}{1.167 - 0.210} = \frac{0.245}{0.957}$$

$$K_2 = \frac{1}{0.02} \times \frac{1.412 - 1.010}{1.010 - 0.210} = \frac{0.402}{0.800} \times \frac{1}{0.02} = 25.12$$

$$K_4 = \frac{1}{0.04} \times \frac{1.412 - 0.817}{0.817 - 0.210} = \frac{1}{0.04} \times \frac{0.592}{0.607} = 24.4$$

$$K_6 = \frac{1}{0.06} \times \frac{1.412 - 0.781}{0.701 - 0.210} = \frac{1}{0.06} \times \frac{0.711}{0.491} = 24.1$$

$$K_8 = \frac{1}{0.08} \times \frac{1.412 - 0.551}{0.551 - 0.210} = \frac{1}{0.08} \times \frac{0.799}{0.403} = 24$$

$$K_{10} = \frac{1}{0.10} \times \frac{1.412 - 0.551}{0.551 - 0.210} = \frac{1}{0.10} \times \frac{0.861}{0.341} = 25.21$$


$$K_{13} = \frac{1}{0.13} \times \frac{1.412 - 0.498}{0.498 - 0.210} = \frac{1}{0.13} \times \frac{0.914}{0.288} = 24.41$$

$$K_{15} = \frac{1}{0.15} \times \frac{1.412 - 0.461}{0.461 - 0.210} = \frac{1}{0.15} \times \frac{0.951}{0.251} = 25.25$$

$$K_{17} = \frac{1}{0.17} \times \frac{1.412 - 0.443}{0.443 - 0.210} = \frac{1}{0.17} \times \frac{0.969}{0.233} = 24.46$$

$$K_{20} = \frac{1}{0.2} \times \frac{1.412 - 0.417}{0.417 - 0.210} = \frac{1}{0.2} \times \frac{0.995}{0.207} = 24.6$$

$$K_{22} = \frac{1}{0.22} \times \frac{1.412 - 0.401}{0.401 - 0.210} = \frac{1}{0.22} \times \frac{1.011}{0.191} = 24.0$$



(iv.) C_t value is taken from the R-solution after one minute two minute or gap of minute till 30 minute by the help stop watch and conductometer.

(v.) to taken by the C-value.

observation table:-

Time	Conductance millimhos	$C_0 - C_t$ millimhos	$C_t - C_{\infty}$	$\frac{C_0 - C_t}{C_t - C_{\infty}}$	K
101	1.412				
1.	1.617	0.245	0.957	0.256	25.60
2.	1.010	0.402	0.800	0.502	25.12
4.	0.817	0.592	0.607	0.980	24.20
6.	0.701	0.711	0.491	0.448	24.13
8.	0.613	0.799	0.403	1.982	24.78
10.	0.551	0.861	0.341	2.521	25.24
13.	0.498	0.914	0.288	3.172	24.41
15.	0.461	0.951	0.251	3.791	25.25
17.	0.443	0.969	0.233	4.163	24.46
20.	0.417	0.995	0.207	4.84	24.03
22.	0.407	1.011	0.191	5.291	25.31

Teacher's Signature : _____

$$K_{25} = \frac{1}{0.25} \times \frac{1.412 - 0.374}{0.374 - 0.210} = \frac{1}{0.25} \times \frac{1.038}{0.164} = 25.3$$

$$K_{27} = \frac{1}{0.27} \times \frac{1.412 - 0.369}{0.369 - 0.210} = \frac{1}{0.27} \times \frac{1.043}{0.159} = 24$$

K_{avg}

$$25.6 + 25.12 + 24.5 + 24.13 + 24.78 + 25.24 + 24.41 + 25.25 \\ + 24.46 + 24.03 + 24.05 + 25.31 + 24.29 + 25.6 \\ = 346.87$$

$$K_{avg} = \frac{346.87}{14} = 24.77 \text{ mol}^{-1} \text{ min}^{-1}$$

$$K \text{ from graph} = \frac{1}{0.01} \times \frac{C_0 - C_t}{C_t - C_{\infty}}$$

$$= \frac{2.3}{9.5 \times 0.1} = 24.21 \text{ mol}^{-1} \text{ min}^{-1}$$

Date

Expt. No.

Page No.

25.	0.374	1.038	0.164	6.333	25.31
27.	0.369	1.048	0.159	6.562	24.29
30.	0.348	1.064	0.138	7.711	25.70
∞	0.210				

Result:-

The rate constant for saponification reaction of ethyl acetate is obtained 24.77 mol⁻¹ min and by graph 24.21 mol⁻¹ min⁻¹

Suresh
23/4/18

Object :-

Determine the triangular phase diagram of acetic acid, chloroform, & water.

Principle :-

Chloroform & water are the practically miscible liquids only one binodal when corresponding to this pair of liquids will be of appear in the phase diagram.

But chloroform is miscible in acetic acid & water mixture of these miscible component are ~~pure~~ pure and to each as the mixture a known amount of third component H₂O is added attentively known amount of two component are mixed and quantities of the third are added until the system to ~~ingentous~~ the quantities of the three component for which this occurs give one point on the binodal curve this complete binodal never is obtained.

Requirements :-

Small glass stoppered bottles
benzite
acetic acid
thoat form.

Teacher's Signature : _____

procedure :-

(1) first mixture of chloroform in acetic acid in stopped bottle with the percentage of 10%, 20%, 30%, 40%, 50%, 60%, 70%, 80% & 90% are prepared.

Now take 20ml CHCl_3 (10%) in a flask and take with H_2O which in m sat H_2 and record the amount of water added when turbidity just appeared.

Above process is done in similar manner the other mixtures.

- calculate the percentage by weight at each component when turbidity just appears.

observation table :-

Amount CHCl_3 (ml)	Amount CH_3COOH	Volume of H_2	Total	H_2O	CHCl_3	CH_3COOH
1	9	1.8	11.8	15.25	8.471	76.27
2	8	2.3	12.3	18.69	16.26	65.04
3	7	2.7	12.7	21.25	23.62	55.11
4	6	3.0	13.0	23.07	30.76	46.15
5	5	3.5	13.0	25.92	37.06	37.03
6	4	3.1	13.5	23.66	45.88	30.54
7	3	2.8	12.8	21.87	54.08	23.43

Teacher's Signature : _____

Calculation:-

(1) for 10% CHCl₃ -

(a) % H₂O = $\frac{1.8}{11.8} \times 100 = 15.25$

(b) % CHCl₃ = $\frac{1}{11.8} \times 100 = 8.471$

(c) % CH₃COOH = $\frac{9}{11.8} \times 100 = 76.27$

(2) for 20% CHCl₃ -

(a) % H₂O = $\frac{2.3}{12.3} \times 100 = 18.69$

(b) % CHCl₃ = $\frac{2}{12.3} \times 100 = 16.26$

(c) % CH₃COOH = $\frac{8}{12.3} \times 100 = 65.04$

(3) for 30% CHCl₃ -

(a) % H₂O = $\frac{2.7}{12.7} \times 100 = 21.25$

(b) % CHCl₃ = $\frac{3}{12.7} \times 100 = 23.62$

(c) % CH₃COOH = $\frac{7}{12.7} \times 100 = 55.11$

(4) for 40% CHCl₃ -

(a) % H₂O = $\frac{3}{13.0} \times 100 = 23.07$

(b) % CHCl₃ = $\frac{4}{13} \times 100 = 30.76$

(c) % CH₃COOH = $\frac{6}{13} \times 100 = 46.15$

(5) for 50% CHCl₃ :-

(a) % H₂O = $\frac{3.5}{13.5} \times 100 = 29$

(b) % CHCl₃ = $\frac{5}{13.5} \times 100 = 37$

(c) % CH₃COOH = $\frac{5}{13.5} \times 100 = 37$

(6) for 60% CHCl₃ :-

(a) % H₂O = $\frac{3.1}{13.1} \times 100 = 23$

(b) % CHCl₃ = $\frac{6}{13.1} \times 100 = 45$

(c) % CH₃COOH = $\frac{4}{13.1} \times 100 = 30$

(7) for 70% CHCl₃ :-

(a) % H₂O = $\frac{2.8}{12.8} \times 100 = 21.8$

(b) % CHCl₃ = $\frac{7}{12.8} \times 100 = 54$

(c) % CH₃COOH = $\frac{3}{12.8} \times 100 = 23$

(8) for 80% CHCl₃ :-

(a) % H₂O = $\frac{2.2}{12.2} \times 100 = 18$

(b) % CHCl₃ = $\frac{8}{12.2} \times 100 = 65$

(c) % CH₃COOH = $\frac{2}{12.2} \times 100 = 16$

(9) for 90% CHCl₃ :-

(a) % H₂O = $\frac{1.7}{11.7} \times 100 = 14.5$

(b) % CHCl₃ = $\frac{9}{11.7} \times 100 = 76$

(c) % CH₃COOH = $\frac{1}{11.7} \times 100 = 8$

Expt. No.

80	8	2	2.2	12.2	18.03	65.57	16.39
90	9	1	1.7	11.7	14.52	76.92	8.54

Result:-

Triangular phase diagram for acetic acid chloroform and water is formed an bimodal area is formed.



Teacher's Signature : _____

Object:-

To determine the strength of acids from given acidic mixture of two acids using NaOH solution by conductometrically.

Principle:-

When a mixture containing acetic acid and hydrochloric acid is titrated with an alkalic strong acid (HCl) will be neutralized first the neutralization of weak acid conductance titration curve will be marked by two points the first and second beaker respectively then

**Requirement:-**

conductivity meter microbottle pipette
breaker.

Chemicals:-

$\frac{N}{10}$ NaOH standard solution and unknown

CH_3COOH unknown HCl

Teacher's Signature : _____

procedure:-

- (i) burette is fill with the standard $\frac{N}{10}$ NaOH.
- (ii.) 10ml HCl and 20ml CH_3COOH is taken in 100ml beaker.
- (iii.) Now set the conductometer on check.
- (iv.) conductometer is set on cell constant.
- (v.) after the conductivity of acid to be checked.
- (vi.) Now 0.5ml NaOH from the burette is added to this solution and check the conductivity this process is repeated till retification result is come and each successive addition of NaOH.

observation:-

Strength of standard NaOH solution $\frac{N}{10}$
Volume of acid 10ml HCl, 20ml CH_3COOH

Volume of NaOH for H, solution $V_1 = 5\text{ml}$

Volume of NaOH for CH_3COOH solution = $V_2 = 1$

$$= 9.9 - 5 = 4.9\text{ml}$$

Calculation:- for HCl:-

$$N_1 V_1 = N_2 V_2$$

$$N_1 = \frac{N_2 V_2}{V_1}$$

$$= \frac{1}{10} \times \frac{5}{30}$$

$$= 0.0167 N$$

$N_1 =$ Normality of HCl = ?
 $V_1 =$ Volume of HCl = 30 ml
 $N_2 =$ Normality of NaOH = $\frac{1}{10} N$
 $V_2 =$ Volume of NaOH = 5 ml

Concentration of HCl in given mixture

$$= N_1 \times \text{equivalent weight}$$

$$= 0.0167 \times 36.5$$

$$= 0.06095 \text{ gm/l}$$

for CH_3COOH :-

$$N_3 V_3 = N_4 V_4$$

$$N_3 = \frac{N_4 V_4}{V_3}$$

$$= \frac{1}{10} \times \frac{4.9}{30} N$$

$$= 0.1633 N$$

concentration of CH_3COOH in given mixture

$$= 0.1633 \times 60$$

$$= 0.9800 \text{ gm/l}$$

Expt. No.

Date
 Page No.

observation table:-

Sr. No.	Volume of alkali (NaOH) addition (ml)	conductance of the solution (milimho)
1.	0.0	2.12
2.	0.5	2.01
3.	1.0	1.821
4.	1.5	1.612
5.	2.0	1.725
6.	2.5	1.592
7.	3.0	1.426
8.	3.5	1.351
9.	4.0	1.296
10.	5.0	1.027
11.	6.0	1.101
12.	7.0	1.151
13.	8.0	1.223
14.	9.0	1.251
15.	10.0	1.313
16.	10.5	1.648
17.	11.0	2.002
18.	11.5	2.398
19.	12.0	2.653
20.	12.5	3.025
21.	13.0	3.412

Teacher's Signature :

Date

Expt. No.

Page No.

Result:-

The strenght of hydrochloric acid is 0.6095
gm/litter and acetic acid is 0.9800 gm/litter
in given acidic mixture.

object :-

calibration of volumetric glass ware.

Introduction :-

calibration is the process by which a stated measure such as the volume of container is checked for accuracy.

- In general, measurement of mass can be determine more precisely and accurately the measurement of volume.
- There for the mass of liquid contained or dispersed by the glass ware will be measured and the corresponding volume calculated using the density of liquid.
- How ever a relatively small change in temp. causes a change in the liquid's volume and thus it's density.

Theory :-

The result of volumetric analysis are based and accrued measurement of volumes of the solution the react with one another. To achieve accurate one mixture employ accurate measuring tools. burettes, pipettes range and variety are used in normal volumetric analysis.

- Therefore pipettes, burettes and measuring cylinder shall be properly calibrated to achieve accurate result.

- calibration of pipettes, burettes and measuring cylinder be done by measuring weight of water delivered. particular glassware.
- Then with the density of water of the volume of 1gm of water at a measured temp. the correct volume is calculated.
- If w is the weight water delivered then

$$\text{Volume water delivered} = w \times v (\text{mL})$$

where v is the volume of 1gm of water in mL at $t^\circ\text{C}$.

Basic principle of glassware calibration:-

- small amount of error in calibration inevitable.
- Tolerance -
 - How much error is allowed in the calibration of a volume measuring item.
- most accurately calibration glassware termed volumetric
- Has narrowest tolerance allowed.
- Two effect of temp. to consider in high accurate measurement.
 - Expand or contract with change in temp.
 - Affects accuracy.
- Volume of the solution change as the temp. changes.
 - Devices calibrated at 20°C with water.

Teacher's Signature : _____

calibration of volumetric flask :-

- weight accurately a thoroughly cleaned and dried flask on a robust balanced.
- Fill the flask with ore free distilled water, so that the lower edge of the meniscus stands at the fixed mark of the neck. Remove any drop of water above the mark by a piece of filter paper.
- dry the outer surface and weight the flask again. After having determined the weight of water contained in the flask upto the mark, obtain the true volume of the vessel from the table → In case the error is considered etch a new ring on the neck.

Temperature °C	Apparent wt of 1cm ³ of water (g)	vol. corresponding to an apparent wt of 1g water (cm ³)	Temp. (°C)	Apparent wt of cm ³ of water (g)	volume on apparent wt of 1g water (cm ³)
10	0.9986	1.0013	18	0.9976	1.0024
11	0.9985	1.0014	19	0.9974	1.0026
12	0.9984	1.0015	20	0.9972	1.0028
13	0.9983	1.0017	21	0.9970	1.0030
14	0.9982	1.0018	22	0.9967	1.0033
15	0.9981	1.0019	23	0.9965	1.0035
16	0.9979	1.0021	24	0.9963	1.0037
17	0.9977	1.0023	25	0.9960	1.0040

Surendra
23/11/18

Teacher's Signature : _____